# Product datasheet Anti-CD45/PTPRC Antibody Catalog Number: A00555-4

BOSTER antibody and ELISA experts

**BOSTER BIOLOGICAL TECHNOLOGY** 

Special NO.1, International Enterprise Center, 2nd Guanshan Road, Wuhan, China

Web: www.boster.com.cn Phone: +86 027-67845390 Fax: +86 027-67845390 Email: boster@boster.com.cn

<b>Basic Inform</b>	nation	
Product Name	Anti-CD45/PTPRC Antibody	
Gene Name	PTPRC	
Source	Rabbit	
Isotype	IgG	
Species Reactivity	mouse, rat	
Tested Application	WB, IHC, ICC/IF, ELISA	
Contents	500 ug/ml antibody with PBS ,0.02% NaN3 , 1 mg BSA and 50% glycerol.	
Immunogen	E.coli-derived mouse CD45/Ptprc recombinant protein (Position: N209-E1269).	
concentration	500 ug/ml	
Purification	Immunogen affinity purified.	
Observed MW	180-250KD	
Dilution Ratios	Western blot (WB): Immunohistochemistry in paraffin section (IHC): Immunocytochemistry/Immunofluorescence (ICC/IF): ELISA: (Boiling the paraffin sections in 10mM citrate buffer,pH6.0,mins is required for the staining of formalin/paraffin section must be determined by end user.	

#### **Storage**

12 months from date of receipt,  $-20^{\circ}$ C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.

## **Background Information**

CD45(Cluster of Differentiation 45), also known as PTPRC, LCA or CD45R, is an enzyme that, in humans, is encoded by the PTPRC gene. CD45 is a member of the protein tyrosine phosphatase (PTP) family. CD45 is a major high molecular mass leukocyte cell surface molecule which is also an integral membrane protein tyrosine phosphatase. The cytogenetic location of CD45 is 1q31.3-q32.1. CD45 is especially a prototype for transmembrane protein-tyrosine phosphatase (PTP). Targeted disruption of the CD45 gene leads to enhanced cytokine and interferon receptor-mediated activation of JAKs and STAT proteins. In vitro, CD45 directly dephosphorylates and binds to JAKs. Functionally, CD45 negatively regulates interleukin-3-mediated cellular proliferation, erythropoietin-dependent hematopoiesis, and antiviral responses in vitro and in vivo. CD45 has been best studied in T cells, where it determines T cell receptor signaling thresholds. CD45 is moved into or out of the immunological synapse (IS) membrane microdomain depending on the relative influence of interaction with the extracellular galectin lattice or the intracellular actin cytoskeleton. Galectin interaction can be finetuned by varying usage of the heavily Oglycosylated spliced regions and sialylation of Nlinked carbohydrates.

### **Selected Validation Data**

#### **Product datasheet**

#### **Anti-CD45/PTPRC Antibody**

Catalog Number: A00555-4



**BOSTER BIOLOGICAL TECHNOLOGY** 

Special NO.1, International Enterprise Center, 2nd Guanshan Road, Wuhan, China

Web: www.boster.com.cn Phone: +86 027-67845390 Fax: +86 027-67845390 Email: boster@boster.com.cn




Figure 1. Western blot analysis of anti- CD45/PTPRC Antibody (A00555-4). The sample well of each lane was loaded with 50ug of sample under reducing conditions.

Lane 1: mouse spleen tissue lysates,

Lane 2: mouse thymus tissue lysates.

Use rabbit anti- PTPRC 1:1000, probed with a goat anti-rabbit IgG-HRP secondary antibody. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002). A specific band was detected for PTPRC at approximately 220KD. The expected band size for PTPRC is at 145KD.