Product datasheet

Anti-LAMC1/LAMC2/LAMC3 Antibody

Catalog Number: BA1761-1



BOSTER BIOLOGICAL TECHNOLOGY

Special NO.1, International Enterprise Center, 2nd Guanshan Road, Wuhan, China

Web: www.boster.com.cn Phone: +86 027-67845390 Fax: +86 027-67845390 Email: boster@boster.com.cn

Basic Information	
Product Name	Anti-LAMC1/LAMC2/LAMC3 Antibody
Gene Name	Lamc1/Lamc2/Lamc3
Source	Rabbit
Isotype	IgG
Species Reactivity	mouse,rat
Tested Application	IHC
Contents	500 ug/ml antibody with PBS $_{\odot}$ 0.02% NaN3 , 1 mg BSA and 50% glycerol.
Immunogen	Human fibroblasts derived.
concentration	500 ug/ml
Purification	Immunogen affinity purified.
Dilution Ratios	Immunohistochemistry in paraffin section (IHC): 1:50-400 (Boiling the paraffin sections in 10mM citrate buffer,pH6.0,or PH8.0 EDTA repair liquid for 20 mins is required for the staining of formalin/paraffin sections.) Optimal working dilutions must be determined by end user.

Storage

12 months from date of receipt, -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.

Background Information

Laminins are major proteins in the basal lamina (one of the layers of the basement membrane), a protein network foundation for most cells and organs. Laminins form independent networks and are associated with type IV collagen networks via entactin, fibronectin, and perlecan. They are important and biologically active parts of the basal lamina, influencing cell differentiation, migration, and adhesion, as well as phenotype and survival. Laminins are trimeric proteins that contain an α -chain, a β -chain, and a γ -chain, found in five, four, and three genetic variants, respectively. Laminins critically contribute to cell attachment and differentiation, cell shape and movement, maintenance of tissue phenotype, and promotion of tissue survival.

Reference

Anti-LAMC1/LAMC2/LAMC3 Antibody被引用在6文献中。

Selected Validation Data

Product datasheet

Anti-LAMC1/LAMC2/LAMC3 Antibody

Catalog Number: BA1761-1



BOSTER BIOLOGICAL TECHNOLOGY

Special NO.1, International Enterprise Center, 2nd Guanshan Road, Wuhan, China

Web: www.boster.com.cn Phone: +86 027-67845390 Fax: +86 027-67845390 Email: boster@boster.com.cn

